

# MOLDCARE® Cushion

*Provides customized patient support*



MOLDCARE® Cushions provide customized head and neck support. Patient comfort is a critical component of patient positioning. The MOLDCARE cushion is composed of a soft fabric that contains expanded polystyrene beads coated in a moisture-cured resin. When activated with water, the cushion becomes rigid, conforming to the contours of the patient's head and neck. It can be used alone or on a Silverman Head Support to create the most comfortable and repeatable head and neck patient specific support.



**RT-4492S**  
Small, 15 cm x 20 cm



**RT-4492**  
20 cm x 25 cm



**RT-4492T**  
Trapezoid, Head &  
Shoulder, 58 cm x 65 cm



**RT-4492C**  
Head & Shoulder S-Frame



**RT-4492U**  
20 cm x 35 cm  
(Shown molded & with Q2 Headrest)



**RT-4492-B**  
Spray Bottle

<b>RT-4492</b>	MOLDCARE® Head Cushion, 20 cm x 25 cm
<b>RT-4492M</b>	MOLDCARE® Head Cushion, 20 cm x 45 cm
<b>RT-4492S</b>	MOLDCARE® Head Cushion, Small, 15 cm x 20 cm
<b>RT-4492T</b>	MOLDCARE® Trapezoid Head & Shoulder Cushion, 58 cm x 65 cm
<b>RT-4492L</b>	MOLDCARE® Cushion, Large, 40 cm x 60 cm
<b>RT-4492XL</b>	MOLDCARE® Cushion, Extra-Large, 60 cm x 60 cm
<b>RT-4492U</b>	MOLDCARE® U Head Cushion, 20 cm x 35 cm
<b>RT-4492C</b>	MOLDCARE® Head and Shoulder S-Frame Cushion
<b>RT-4492-B</b>	Spray Bottle for MOLDCARE®



# Silverman Head Support

*Provides customized patient support*



Silverman head and neck supports are available in a variety of sizes to fit your positioning needs.

- The Silverman Q1 and Q2 minimize overstretching of the thermoplastic mask while simultaneously offering cranial support
- Silverman & MOLDCARE together are lower attenuating than most standard foam headrests
- Q Headrests feature a posterior cranial indexing notch for enhanced reproducibility of MOLDCARE cushions
- Q Headrests extend to the top of the cranium for added support

## Silverman Q Head Supports



**RT-4491-Q1**  
Silverman Q1



**RT-4491-Q2**  
Silverman Q2



**RT-4491-Q3-2G**  
Silverman Q3



**RT-4491-Q4-2G**  
Silverman Q4



**RT-4491R-M\***  
Silverman 'M' Low-Profile  
Head Support



\*Exclusive low-profile design allows the occiput to rest below the top plane of the board. By lowering the patient's head, the mask is stretched a shorter distance during the molding process. Less stretch means a more comfortable mask and an overall better fit for your patient.

## Silverman Head Support



**RT-4491R2-A**  
size A



**RT-4491R2-B**  
size B



**RT-4491R2-C**  
size C



**RT-4491R2-D**  
size D



**RT-4491R2-E**  
size E



**RT-4491R2-F**  
size F

**RT-4491R2** Silverman Head Supports, set of 6, sizes A-F